TECHNO LOGY AND SERVICES CATAL OGUE

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**DUNA'S INNOVATION
& TECHNOLOGY CONCEPT**

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Duna has always believed in technology to improve its production processes. The company has been active in the field of process innovation addressed to design and manufacture customized orthopaedic footwear and orthosis for more than 20 years. Thanks to Duna's partnership with important public and private research centers in challenging projects, the possibility of having a built-in technology to support these processes is now a reality.

“Technology thus becomes the centre of the company's customer added-value services. It is the turning point to improve the ability to satisfy not only the customer's demands in terms of aesthetic rendering but also and especially in terms of walking functionality and therapeutic purposes.”

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ORTOSHOE

INNOVATION IN MANUFACTURING OF CUSTOMIZED ORTHOPAEDIC SHOES

2009 >

SPECIAL SHOES MOVEMENT

GRANT AGREEMENT NO:
 NMP - 2008 - SME -
 2 - R: 229261

2011 >

TCO

TECHNOLOGICAL CUSTOMIZED ORTHOSIS

2013 >

INPRUF

INNOVATIVE PRODUCTION PROCESS FOR MULTIFUNCTIONAL CUSTOMIZED SHOES

2015 >>

In more than 40 years, Duna has constantly increased the range of orthopaedic products thanks to its creativity, innovative behaviour and production quality to grant the design, manufacturing and sales of insole-preset footwear and customized footwear/orthosis for any age and with specific characteristics to treat any foot disease and patient's posture. In partnership with Otto Bock since 1994, a world leader in healthcare, Duna is today national market leader, one of the major European players in the market, and operates directly or through its distributors in more than 12 countries.

In 2008 Duna opened the Innovation & Technology department with the aim of developing technology products and added-value services to enhance the interactive processes between orthopaedics and production centers. This provided customized shoes and orthotics able to satisfy the end-users and to improve their quality of life. The technological products and services developed directly by Duna's Innovation & Technology Dept. are the result of research and development oriented to overcome the limitations of those systems created for serial production, to realize integrated systems to optimize processing performance, to produce modular systems "opened" to other technologies to guarantee the maximum flexibility in each orthopaedic shop. These systems are used in everyday processes where we also play as technological provider and customized product supplier. "A customized device is not a product but a service focused on the patient's needs and we are able to understand your needs because these are our daily needs".

Duna's Innovation & Technology shows three "technology programs" with a product portfolio that can satisfy any need. Consulting, training and after-sales technical support is carried out by highly qualified personnel.

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

TRAINING
SCHULUNG

MATERIALS
MATERIALIEN

TECHNOLOGY FOR CUSTOMIZED FOOTWEAR

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OrthoShoe process is the summary of Duna's technology vision in terms of customization of orthopaedic footwear. The hardware and software tools presented by Duna's Innovation & Technology have been designed and are constantly developed to combine modularity and full integration with the aim of improving the workflow between the orthopedic technician or podiatrist and the production center and to provide the final patient with a higher quality product in a short time. OrthoFootScan2 and FootWizard® ensure a fast and complete control of the digital acquisition process of the foot, the technical configuration of the shoe and data sending to the production center. The CAD sw application Last Designer is intended to be used by people trained for the stylization of orthopaedic last. Orthoshoe is based on a shared approach towards the patient and can be applied completely or partially to fit the needs of each technician and organization.

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


3D DIGITAL ACQUISITION OF FOOT MEASURES

3D GRÖSSEN ERWERB

OrthoFootScan 2

The 3D full foot digitalization is really important because it avoids mistakes in the acquisition of manually taken measurements and allows to make an orthopaedic customized last keeping the same accuracy of a last made using a foot cast. Thanks to a 3D approach, the workflow is faster, avoids health risks to the worker and transport costs to the production center.

The transportable 3D **OrthoFootScan2** proposed by Duna's Innovation & Technology is based on mixed stereo-cameras and blue laser technology, allowing a very fast and high level of accuracy in the digitalization also for the most serious pathological feet. The system is additionally equipped with a heel bearing support to resemble the same heel height the foot will have inside the shoe. The **Foot3D** software allows the user to take the notable girths across the foot in a very friendly way and to personalize the placement of the measurement planes with few clicks. The acquisition in *.stl file format makes the system opened towards any CAD platform while the measures taken can be uploaded automatically into **FootWizard**[®].

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
TECHNICAL FOOTWEAR CONFIGURATION AND ORDER MANAGEMENT TO THE PRODUCTION CENTER

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Technical configuration of the footwear and the completeness of data sent to the production company are important aspects as much as the precision of the acquisition system in order to be sure to avoid non-conformities in the manufacturing of the final shoe.

FootWizard® is a knowledge based software that represents Duna's long-lasting experience in the field of customized orthopaedic footwear and makes it at disposal of the patient's well-being.

FootWizard® Custom totally supports the orthopaedic technician in technical and administrative management of the custom-made footwear. It allows to monitor all the necessary information, to reduce mistakes and wrong interpretation, to give help to technical and aesthetical choices, to speed up the workflow and reduce operative costs. Operating in more than 10 languages, few simple steps are necessary in **FootWizard® Custom** to manage administrative data of the order, to insert foot measurements, fill in the technical order and watch a realistic 3D rendering of the final custom-made footwear.

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
ORTHOPAEDIC CUSTOMIZED LAST DESIGN

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Last Designer

Shoe last is the most important structural element in the manufacturing of footwear. Starting from foot measurement, Last Designer and styling is necessary to satisfy both functional and aesthetical requirements. In case of foot diseases, the last design process is even more important in comparison with respect to the traditional footwear industry, because the last is the structure that also defines the corrections and the orthopaedic performance of the final product.

Software application **Last Designer** is based on CAD Rhinoceros plug-in developed by Duna's Innovation & Technology. It implements automatic algorithms that are the synthesis of the experience acquired by Duna's shoe last craftsmen and makes a CAD last design possible starting from the acquisition of a cloud-point through a 3D scanner, or manually taken measures or directly from a positive foot cast.

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
TECHNOLOGY FOR CUSTOMIZED INSOLES

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
TCO – Technological Customized Orthosis – it is the innovative measure-to-manufacturing production process developed by Duna's Innovation & Technology to innovate all the operating phases that contribute to the manufacturing of orthopaedic custom-made insole, even if it is used in regular fitting or individual shoes. HW e SW tools have been thought to provide the final user with a system capable to adapt to all the possible cases of daily work with efficiency and functionality, reducing to zero the necessity to operate manually on the finished product and optimizing the entire production process (crossing times and production costs). The key element of the whole process is represented by the TCO Insole Designer, a 3D CAD for designing, able to integrate, even during the planning phase, the insole and the orthopaedic footwear that will guest it, maximizing the therapeutic effect and the final quality for both the operator and especially the end user. TCO is therefore a method that the Orthopaedic shop is free to perform independently or through an outsourcing scheme (design and/or production), thus defining the more suitable production process according to its necessities. The choice to operate on opened file formats has been made to allow the users to decide which and how many tools of the TCO process they might introduce in their shop.

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Foot digitalization of the patient is the starting step to join a CAD-CAM productive workflow and it must be considered equal to foot cast styling or phenolic foam footprint capturing. Duna's Innovation & Technology offers different possible technological solutions to the end-user to reach this task efficiently and reliably.

3D DigiPlan scanning machine allows to obtain the bottom of the foot on a transparent flat surface, ensuring a high measurement accuracy and the additional advantage represented by the portability of the whole system. **OrthoFootScan2** machine adds the advantage represented by the possibility to place the foot at the correct heel height, avoiding to deform the scanning inside the CAD and thus to make dimensional mistakes on the final product. At the same time, this acquisition can be used for the designing both the insole and the shoe, integrating the two products. The **TCO Scan** system also provides a full contact with the foot and enables the foot to be detected as it were placed inside the final footwear (also including possible corrections).

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
TECHNICAL CONFIGURATION AND ORDER MANAGEMENT TO THE PRODUCTION CENTER

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The complete management of the relevant information is very important to ensure successful design and construction of the insole at the production center both for internal and external production solutions.


FootWizard® is a knowledge based software that introduces Duna's long-lasting experience in the field of customized orthosis for personalised orthopaedic footwear, making it available for the patient's well-being.

FootWizard® Insole allows to manage the patient/orders database and, in few simple steps, to convey the technical data to the production center including the foot scanning or the project file related to the insole with the additional advantages of the traceability and completeness of information.

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Insoles addressed to people affected from foot diseases are medical devices and, together with the shoes where they are placed, have to contribute to restore gait efficiency, reduce energy consumption and finally produce a biomechanical and therapeutic effect for the patient. Their structure has to satisfy the patient's needs and not to be constrained to predefined structural limits.

TCO Insole Designer software application developed by Duna's Innovation & Technology allows to automate the design of custom-made insoles, considering the final shoe used. The internal database allows to reach the final result in few simple clicks, operating on a set of constantly updated templates. In complex cases it is possible to adjust the insole starting from the patient's foot, without any restriction or limit to the structures and achievable functionality. Developed as plug-in of CAD Rhinoceros, it carries out a complete operating bundled together with Last Designer software and CAM software for the tool path generation.

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INSOLE PRODUCTION

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MILLING
MACHINES
PANTOGRAPHEN

The automatic production of a customized insole requires the use of CNC controlled milling machines, designed and optimized for working on expanded materials. The passage from a digital CAD model to a physical realization requires the use of a specific software (CAM) to generate the path for the tool equipping the numerical control machine. The possibility to set the main configuration parameters of the machine in a flexible way, such as feed-rate or spindle spin-rate, in addition to the experience in the use of the right materials and tools, directly affects production time and finishing quality.

Duna's Innovation & Technology proposes different solutions to meet each customer's operational requirements. The available systems are extremely strong and designed to make any kind of product in a wide range of possible materials and geometries. Each system can be customized with a long list of options to satisfy the needs of both a small and large production. Through our internal staff, we are able to provide pre-sale technical consulting services, after-sale assistance, supply materials and replace consumables.

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ADVANTAGES AND SERVICES OF DUNA'S TECHNOLOGY

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DUNA
INNOVATION
& TECHNOLOGY

Why should we use DUNA's technology in footwear and insole customization?

All the technological solutions are designed to be fully integrated one with another but at the same time 'open' to be used in already partially digitized processes

All hardware and software systems proposed are the result of research and development activities carried out with important national and international research centers. They have been integrated and tested by Duna's staff and are applied in everyday production of orthopedic customized footwear and insole.

Thanks to its highly qualified and specialised staff, Duna's technological innovation, is able to provide:

TECHNICAL CONSULTING during the pre-sale phase to identify solutions that ensure the best performance and the best feedback on investments.

TRAINING at the customer's or at Duna's headquarters with regard to the use of hardware and software tools provided, order and configuration processes of customized products, use of materials and the integration capacity coming from the assessment of biomechanical performance using foot and gait diagnostic systems.

AFTER-SALE TECHNICAL ASSISTANCE accomplished directly by Duna's people to ensure constant support and the maximum speed in the case of intervention at the customer's.

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FootWizard®

Foot Wizard® R4 is a unique system to ensure a complete technical configuration and management of the order in the production center for the manufacturing of customized footwear and insoles. Suitable algorithms and a powerful graphical engine allow the virtual mockup of the shoe, presenting a 3D realistic visualization, starting from individually taken measurements. **Foot Wizard® R4** allows detailed definition for the end-product features, respecting the requirements of the single patient for footwear and insoles, and providing the orthopaedic technician or the podologist with an essential instrument for daily work. The system currently operates in 11 languages, implementing a common technical alphabet between the users, who can employ it as a dialogue tool, with an additional advantage of details and graphical explanations for immediate understanding. In case of multi-shop realities, additionally, the architecture of the system allows the synchronization of internal database of each node with the others, thus simplifying data exchange.

Foot Wizard® R4 is a platform that can be equipped with one or more modules, adapting itself to the needs of each customer:

- _ **Custom** for individual personalization
- _ **Modular** for modular shoes personalization
- _ **Insole** for insoles personalization

MINIMUM PC REQUIREMENTS

Intel Core i3 processor
Graphic Card: 1 GB dedicated
2GB RAM
2GB Free space available on hard-disk
Windows 7 32 bit Operating System

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The **FootWizard® Custom** provides users with important functions:

- _ Data archive. It allows the storage of administrative and technical data of the patient's orders, providing a historical and always updated situation.
- _ Acquisition of foot measures. Total and guided measures panel goes beyond the traditional manually compiled measurement sheet, with a set of detailed, pre-collected and organized information, reducing to zero the possibility of compilation errors. In case of **OrthoFootScan2** usage, this operation is included "in a single click" time.
- _ Technical configuration of the shoe. It allows to define any characteristic the user wants about the orthopaedic shoe. Size, materials, structures, details and colors can be defined with an almost endless number of possible combinations (and always constantly updated!).
- _ 3D real time rendering. Thanks to the embedded graphical engine, it enables 3D display of the footwear that will be produced and delivered, in terms of dimensions, colors and structures.
- _ Data management. Digital transmission of the order using email service and in real-time, zero cost and with additional advantage of the traceability. The complete documentation sent is imported and archived digitally, reducing to zero the possibility of human error in the data transfer to the production center (first reason of non-compliances recorded).

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FootWizard® Insole is the module specifically thought and created to manage the production workflow related to the production of the personalised orthotics, thus becoming the ideal tool to manage and organize both the internal production of the Orthopaedic shop and the total or partial outsourcing of the custom-made insole production.

The user interface is simple and very friendly. It is organized in a sequence of boards and allows the user to configure the order in a very fast and error-avoiding way.

The dedicated measurement panel makes the import of acquired files generated by scanning machines very simple and immediate, especially when Duna's **OrthoFootScan2** or **TCO Scan** are used. Alternatively, it is possible to write foot dimensions by means of dedicated tabs.

The section related to production details allows to define any kind of desired parameter for the final insole: from the materials of each component to the additional items which can be configured through specific and interactive options.

Data sending and processing are completely digital, drastically reducing crossing times and data entry errors.

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Ortho FootScan 2

OrthoFootScan2 has been developed to make 3D digitals scan of feet that need customized insole and footwear. Both high scan volume and a mechanical system to place the foot at the same heel height of the final shoe allow to take the measurements in ideal conditions avoiding the use of the cast. The blue laser technology allows to operate in all lighting conditions. **Foot3D** application allows a personalized measurement's plan displacement and gives the possibility to generate files in standard STL format. The system can be comfortably transported in its anti-shock trolley and can be connected via USB to any portable computer which can be positioned inside the trolley. **OrthoFootScan2** interfaces directly with the order configuration tool and virtual mock-up **FootWizard®**.

TECHNICAL FEATURES

Scanning Volume LxHxW = 400x160x140 mm
High Precision Accuracy = +/- 1 mm
Scanning rate: 18 sec
Data processing rate: 30 - 50 sec
Output 3D Format : STL, ASCII
Weight supported: 120 Kg
Dimensions LxHxW: 735 x 416 x 422 mm
Weight: 18 Kg
Power Supply: 100-240Vac 50/60Hz 150W
Connectivity: USB

EQUIPMENT

Foot Scanner
Software Foot3D with USB Safety Dongle
Stainless steel handrail
Transportability kit with lateral supports and antishock trolley LxHxP=735 x 416 x 422mm
Foot platform for fixed placement (optional)

MINIMUM PC REQUIREMENTS

Intel Core i3 processor, 512 MB Graphic Card dedicated, 2GB RAM, 2GB Free space available on hard-disk, Windows 7 32 bit Operating System

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TCO Scan combines both highly precision mechanics and laser technologies to reach a 3D foot plantar digitalization in ideal biomechanical conditions with excellent performances:

- _ Possibility to modify the profile of the foot support surface simulating the loading conditions of the foot inside the shoe
- _ Possibility to simulate the presence of pronation and supination planes to reach a foot scan in adjusting conditions
- _ Fast data capture and processing of the foot bottom surface
- _ Dedicated and expandable database of Insole profile
- _ Interface with **Foot Wizard® R4**
- _ User-friendly Graphical User Interface with editing and uploading of set profile
- _ Compatibility with all CAD/CAM Systems

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TECHNICAL FEATURES

Scanning Area LxW = 322x160 mm
Max profile height setting: 80 mm
UltraHigh precision system: +/-30 micron
Scanning rate (with profile setted): 20 sec
Output 3D Format: points cloud (ASCII)
Weight supported: 120 Kg
Dimensions LxHxW: 570 x 490 x 380 mm
Power Supply: 220V 50 Hz 400W
Connectivity: USB

EQUIPMENT

TCO scanner
TCO SW interface with USB Safety Dongle
System Cabling
Plattform with aluminium handrail and extension for chair displacement LxW = 1100 x 950 (1250) mm (optional)

MINIMUM PC REQUIREMENTS

Intel Core i5 processor, 4GB RAM, 2GB Free space available on hard-disk, Windows 7 32 bit Operating System

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
3D DigiPlan

3D Digiplan is a scanning system developed to keep the digitalization of the foot bottom surface. It guarantees a very friendly use and portability. Able to operate in any lighting condition, the system makes a projection of the foot on a flat transparent surface up 60 mm., thus acquiring all the morphological important structures for the insole design. 3D Digiplan has a routine addressed to the scanning of foam boxes, optimizing the result and avoiding transport costs of the materials. STL file output generated by the system owns the ideal features to be used inside 3D CAD design for insole manufacturing.

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TECHNICAL FEATURES

Scanning Volume LxHxW = 400x60x140 mm
Accuracy = +/- 1 mm
Scanning rate: 18 sec
External Dimensions LxHxW: 725x136x420 mm
Weight: 10 Kg
Power Supply: 100-240Vac 50/60Hz 150W
Connectivity: USB
Equipment: Foot ScannerData acquisition and machine control software

MINIMUM PC REQUIREMENTS

Intel Core i3 processor, 2GB RAM, 2GB Free space available on hard-disk, Windows 7 32 bit Operating System

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Last Designer

Last Designer comes from Duna's necessity to develop a sw application able to optimize internal processes oriented to design customized shoes. This sw application, conceived as a Rhinoceros CAD plugin, gathers Duna's know-how in the field of last creation, optimizing the internal designing processes starting from cloud-points originated from 3D foot digitalization, from manually taken measurements or from plaster cast.

Sw application enables an operator having both CAD and last know-how to design the final last in an easy and fast way using dedicated macros. A continuously updated database gives the possibility to design a customized last easily and using few commands thanks to the lasts designed by Duna starting from the study of specific lasts for foot pathologies.

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
MINIMUM PC REQUIREMENTS

Intel Core i3 processor
RAM: 2 GB
Free space available on hard-disk: 5 GB
Graphic Card: 1 Gb dedicated
Operating System: Windows 7 32 bit
Rhinoceros 4.0 license

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TCO Insole Designer is a tool developed to support a 3D CAD insole design and to allow the end user to create a truly customized object according to the patient's needs, without imposing dimensional or structural limitations of any kind. **TCO Insole Designer** offers the possibility to take into consideration the footwear in which the foot has to be inserted during the designing phase. Whether the shoe is to be a regular fitting or customized one, this tool prevents the user from having to make changes and manual adjustments to the product after CNC milling. In the case the user wants to develop a custom-made insole for insole-preset or modular footwear, it is possible to work according to a streamlined workflow that allows to get the finished object in few clicks, while keeping the ideal characteristics of the orthotic.

The insole can be developed starting from a 3D scan of the patient's foot acquired by the scanning machines proposed by Duna or any other system that generates a standard STL file format and considering static footprints or plantar pressure maps. Developed as a plug-in of Rhinoceros CAD, **TCO Insole Designer** creates a single integrated working environment together with **Last Designer** and CAM software for the generation of tool paths.

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MINIMUM PC REQUIREMENTS

Intel core i3 Processor
RAM: 2 GB
Free space available on hard-disk: 5 GB
Graphic Card: 1 GB dedicated
Operating System: Windows 7 32 bit
Rhinoceros 5.0 license SR9

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MILLING MACHINES Pantographen

CNC 3 axes milling machines, specially designed and built for milling orthotics. The control system can handle CAM toolpaths in ISO G-code format (Open System). Various models are available with different working table areas and working heights under the milling head to cover all production needs starting by clamping minimum of 2 pairs. The basic system can be completed with options to customize the machine to meet specific production requirements (dust extraction of the particle systems, probes for the automatic zero setting, control PAD to speed up the execution of repeated commands, Control Panel Multifunction board machine, automatic tool change system, vacuum clamping system etc.).

TECHNICAL FEATURES

Aluminum frame
Work surface: 700x300mm minimum
Z displacement: 80mm minimum
Control system: MAC3 (mod. S2) / TexComputer (mod. S4 or higher)
Linear guides with recirculating ball 15mm
Toothed belt drive in polyurethane and steel wires with inextensible HTD3M step
Spindle rotation speed: 18,000RPM minimal
Fast forward speed: 15 m/min
Resolution: 0.01 mm
Control electronics: CeNeCe C7 microstepping
Stepper motors to high torque 24HB90 and 30 kg/cm (mod. S2)
Brushless motors 1.4 hp (mod. S4 or higher)

For any clarification please contact our technical staff. We who will help you with the choice of the model and the most suitable options suitable to your needs. Duna is authorized reseller of MecSoft Corporation product RhinoCAM and VisualMill.

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


P-WALK is the baropodometric measurement system which allows you to measure and analyse the plantar pressure and the spatio-temporal gait parameters in a simple and efficient way. It consists of a set of sensor-provided platforms, joined together to form a walkway of personalised length (from 50cm to 4m). **P-WALK** evaluates the static and dynamic support, highlighting pressure overload, postural rotations and dysmetria. In addition to the analysis of pressure, **P-STUDIO** software allows you to:

- _ Study the spatio-temporal gait parameters for a functional assessment of the patient under examination, highlighting the therapeutic effects and moving skills;
- _ Make stabilometric analysis to assess the degree of balance of the subject through the study of the Center of Pressure;
- _ Make comparisons with the population's normal classes highlighting the differences of a subject with respect to the average value.

It is also possible to complete the system with cameras (max 2) for the video recording of the movements and subsequent analysis.

TECHNICAL FEATURES

Dimensions: 675 x 540 x 5 mm
Active area: 480 x 480 mm
Sensor's dimensions: 10 x 10 mm
Sensors number: 2304
Acquisition technology: resistive
Acquisition frequency: 100 Hz
(50 Hz with multiple platforms)
Measure span: 30 - 400 kPa
Connectivity: USB 2.0

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
G-WALK is the ideal technical solution for a quick and objective evaluation of walking disturbances.

It is made of a wireless wearable inertial sensor that allows to measure quantitatively the body kinematic in its principal components: spatio-temporal, symmetry, propulsion and pelvic movement. It is intuitive and easy to use. Tests are extremely quick to be performed and the report generated is completely automatic with an immediate comparison with normative bands.

In the orthopedic field **G-WALK** gives the most important information to plan the rehabilitation after injuries or surgeries, to maximize the movement recovery and to define patient moving capabilities. Also in the prosthetic field and when customizing an orthotic device, **G-G-WALK** shows a wide range of interest allowing O&P specialists to accomplish their work according to the patients' characteristics.

Different modules and test protocols can be added to basic **G-WALK** sensor, increasing the testing capabilities also considering the testing environment.

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© **COPYRIGHT DUNA**
ALL RIGHTS RESERVED

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GRAPHIC DESIGN
GRAPHISCHES DESIGN ADV
COMBINATORIA CREATIVA

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autorizzazione di Duna Srl.
Genehmigung von Duna Srl