

PERSONALIZE OUR SHOES ACCORDING TO YOUR NEEDS.

¡ PERSONALIZA TUS CALZADOS CON TUS NECESIDADES !


Thanks to our innovative tools creating your own personalized pair of shoes is now even easier. Choose the model your prefer from our database, mark the adjustments needed and in few steps your shoes are ready according to the functionality and therapeutic effect required.

Gracias a nuestros nuevos instrumentos está aún más fácil realizar tus zapatos. Elige el modelo que más prefieres en nuestra base de datos, marca las modificaciones necesarias y en pocos pasos tus zapatos estarán listos según la funcionalidad y el efecto terapéutico necesario.



MODEL SELECTION

ELIGE MODELO



COLOR SELECTION

ELIGE COLOR

STIFFENERS ADJUSTMENTS

AJUSTA REFUERZOS

LAST AND WIDTH SELECTION

ELIGE ANCHO Y FORMA

AVAILABLE WIDTHS
8 / 10 / 11 / 12 / 13 / 15

ADJUST THE LAST IN THE GIVEN POINTS

(A. INSTEP / B. BALL / C. TOE HEIGHT)

AJUSTA LA FORMA EN ESTOS
PUNTOS


(A. EMPEINE / B. TERCIO ANTERIOR DEL PIE /
C. PUNTA)

OUTSOLE SELECTION

ELIGE LA SUELA

SELECTION OF DIABETIC LAST MANUFACTURING

ELIGE LA FORMA PARA
ZAPATOS DIABÉTICOS



SHOE DESIGN AFTER ADJUSTMENTS

MODELO DE ZAPATO DESPUÉS
ADAPTACIONES



FINISHED SHOE

ZAPATO FINAL

duna®

ORTHOPAEDIC SHOES INDUSTRY

duna.it